

Antibodies Product Datasheet

Product Name:	Di-Ras1 Antibody
Catalog Number:	C16016
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DIRA1. AA range:149-198
Concentration:	1 mg/ml
Gene Name:	DIRAS1
Protein Name:	GTP-binding protein Di-Ras1
Other Name:	DIRAS1; GBTS1; RIG; GTP-binding protein Di-Ras1; Distinct subgroup of the Ras family member 1; Ras-related inhibitor of cell growth; Rig; Small GTP-binding tumor suppressor 1
Human Swiss Prot No.:	O95057
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Di-Ras1 Polyclonal Antibody detects endogenous levels of Di-Ras1 protein.
Storage and Stability:	-20°C/1 year