

Antibodies Product Datasheet

Product Name:	GUF1 Antibody
Catalog Number:	C16021
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GUF1. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	GUF1
Protein Name:	Translation factor GUF1 mitochondrial
Other Name:	GUF1; Translation factor GUF1; mitochondrial; Elongation factor 4 homolog; EF-4; GTPase GUF1; Ribosomal back-translocase
Human Swiss Prot No.:	Q8N442
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GUF1 Polyclonal Antibody detects endogenous levels of GUF1 protein.
Storage and Stability:	-20°C/1 year