

Antibodies Product Datasheet

Product Name:	Gy 5 Antibody
Catalog Number:	C16046
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GNG5. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	GNG5
Protein Name:	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5
Other Name:	GNG5; GNGT5; Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5
Human Swiss Prot No.:	P63218
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gy 5 Polyclonal Antibody detects endogenous levels of Gy 5 protein.
Storage and Stability:	-20°C/1 year