

Antibodies Product Datasheet

Product Name:	Gα i-3 Antibody
Catalog Number:	C16054
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GNAI3. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	GNAI3
Protein Name:	Guanine nucleotide-binding protein G(k) subunit alpha
Other Name:	GNAI3; Guanine nucleotide-binding protein G(k) subunit alpha; G(i) alpha-3
Human Swiss Prot No.:	P08754
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gα i-3 Polyclonal Antibody detects endogenous levels of Gα i-3 protein.
Storage and Stability:	-20°C/1 year