

Antibodies Product Datasheet

Product Name:	Gα t2 Antibody
Catalog Number:	C16063
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GNAT2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GNAT2
Protein Name:	Guanine nucleotide-binding protein G(t) subunit alpha-2
Other Name:	GNAT2; GNATC; Guanine nucleotide-binding protein G(t) subunit alpha-2; Transducin alpha-2 chain
Human Swiss Prot No.:	P19087
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gα t2 Polyclonal Antibody detects endogenous levels of Gα t2 protein.
Storage and Stability:	-20°C/1 year