

Antibodies Product Datasheet

Product Name:	GBP3 Antibody
Catalog Number:	C16080
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GBP3. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	GBP3
Protein Name:	Guanylate-binding protein 3
Other Name:	GBP3; Guanylate-binding protein 3; GTP-binding protein 3; GBP-3; Guanine nucleotide-binding protein 3
Human Swiss Prot No.:	Q9H0R5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GBP3 Polyclonal Antibody detects endogenous levels of GBP3 protein.
Storage and Stability:	-20°C/1 year