

Antibodies Product Datasheet

Product Name:	GBP5 Antibody
Catalog Number:	C16082
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human GBP5.
Concentration:	1 mg/ml
Gene Name:	GBP5
Protein Name:	Guanylate-binding protein 5
Other Name:	GBP5; Guanylate-binding protein 5; GBP-TA antigen; GTP-binding protein 5; GBP-5; Guanine nucleotide-binding protein 5
Human Swiss Prot No.:	Q96PP8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GBP5 Polyclonal Antibody detects endogenous levels of GBP5 protein.
Storage and Stability:	-20°C/1 year