

Antibodies Product Datasheet

Product Name:	HGFA Antibody
Catalog Number:	C16101
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human HGFA.
Concentration:	1 mg/ml
Gene Name:	HGFAC
Protein Name:	Hepatocyte growth factor activator
Other Name:	HGFAC; Hepatocyte growth factor activator; HGF activator; HGFA
Human Swiss Prot No.:	Q04756
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HGFA Polyclonal Antibody detects endogenous levels of HGFA protein.
Storage and Stability:	-20°C/1 year