

Antibodies Product Datasheet

Product Name:	hnRNP DL Antibody
Catalog Number:	C16111
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HNRPD. AA range:241-290
Concentration:	1 mg/ml
Gene Name:	HNRPD
Protein Name:	Heterogeneous nuclear ribonucleoprotein D-like
Other Name:	HNRPD; JKTBP; Heterogeneous nuclear ribonucleoprotein D-like; hnRNP D-like; hnRNP DL; AU-rich element RNA-binding factor; JKT41-binding protein; Protein laAUF1
Human Swiss Prot No.:	O14979
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	hnRNP DL Polyclonal Antibody detects endogenous levels of hnRNP DL protein.
Storage and Stability:	-20°C/1 year