

Antibodies Product Datasheet

Product Name:	Rabr Antibody
Catalog Number:	C16140
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGFG2. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	AGFG2
Protein Name:	Arf-GAP domain and FG repeat-containing protein 2
Other Name:	AGFG2; HRBL; RABR; Arf-GAP domain and FG repeat-containing protein 2; HIV-1 Rev-binding protein-like protein; Rev/Rex activation domain-binding protein related; RAB-R
Human Swiss Prot No.:	O95081
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rabr Polyclonal Antibody detects endogenous levels of Rabr protein.
Storage and Stability:	-20°C/1 year