

Antibodies Product Datasheet

Product Name:	HLA-H Antibody
Catalog Number:	C16165
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HLAH. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	HLA-H
Protein Name:	Putative HLA class I histocompatibility antigen alpha chain H
Other Name:	HLA-H; HLAH; Putative HLA class I histocompatibility antigen; alpha chain H; HLA-12.4; HLA-AR; MHC class I antigen H
Human Swiss Prot No.:	P01893
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HLA-H Polyclonal Antibody detects endogenous levels of HLA-H protein.
Storage and Stability:	-20°C/1 year