

Antibodies Product Datasheet

Product Name:	HLA-DO α Antibody
Catalog Number:	C16218
Reactivity:	Human
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HLA-DOA. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	HLA-DOA
Protein Name:	HLA class II histocompatibility antigen DO alpha chain
Other Name:	HLA-DOA; HLA-DNA; HLA-DZA; HLA class II histocompatibility antigen; DO alpha chain; MHC DN-alpha; MHC DZ alpha; MHC class II antigen DOA
Human Swiss Prot No.:	P06340
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HLA-DO α Polyclonal Antibody detects endogenous levels of HLA-DO α protein.
Storage and Stability:	-20°C/1 year