

Antibodies Product Datasheet

Product Name:	HLA-DO β Antibody
Catalog Number:	C16219
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HLA-DOB. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	HLA-DOB
Protein Name:	HLA class II histocompatibility antigen DO beta chain
Other Name:	HLA-DOB; HLA class II histocompatibility antigen; DO beta chain; MHC class II antigen DOB
Human Swiss Prot No.:	P13765
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HLA-DO β Polyclonal Antibody detects endogenous levels of HLA-DO β protein.
Storage and Stability:	-20°C/1 year