

## Antibodies Product Datasheet

|                               |                                                                                                                                                                     |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HLA-DP $\alpha$ 1 Antibody                                                                                                                                          |
| <b>Catalog Number:</b>        | C16220                                                                                                                                                              |
| <b>Reactivity:</b>            | Human                                                                                                                                                               |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                           |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.                                                                                 |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                    |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                               |
| <b>Immunogen:</b>             | Synthesized peptide derived from the Internal region of human HLA-DP $\alpha$ 1.                                                                                    |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                             |
| <b>Gene Name:</b>             | HLA-DPA1                                                                                                                                                            |
| <b>Protein Name:</b>          | HLA class II histocompatibility antigen DP alpha 1 chain                                                                                                            |
| <b>Other Name:</b>            | HLA-DPA1; HLA-DP1A; HLASB; HLA class II histocompatibility antigen; DP alpha 1 chain; DP(W3); DP(W4); HLA-SB alpha chain; MHC class II DP3-alpha; MHC class II DPA1 |
| <b>Human Swiss Prot No.:</b>  | P20036                                                                                                                                                              |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                             |
| <b>Specificity:</b>           | HLA-DP $\alpha$ 1 Polyclonal Antibody detects endogenous levels of HLA-DP $\alpha$ 1 protein.                                                                       |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                        |