

Antibodies Product Datasheet

Product Name:	HoxA11/D11 Antibody
Catalog Number:	C16289
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HOXA11/D11. AA range:216-265
Concentration:	1 mg/ml
Gene Name:	HOXA11/HOXD11
Protein Name:	Homeobox protein Hox-A11/D11
Other Name:	HOXA11; HOX11; Homeobox protein Hox-A11; Homeobox protein Hox-1I; HOXD11; HOX4F; Homeobox protein Hox-D11; Homeobox protein Hox-4F
Human Swiss Prot No.:	P31270/P31277
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HoxA11/D11 Polyclonal Antibody detects endogenous levels of HoxA11/D11 protein.
Storage and Stability:	-20°C/1 year