

Antibodies Product Datasheet

Product Name:	HoxA5 Antibody
Catalog Number:	C16294
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human HoxA5.
Concentration:	1 mg/ml
Gene Name:	HOXA5
Protein Name:	Homeobox protein Hox-A5
Other Name:	HOXA5; HOX1C; Homeobox protein Hox-A5; Homeobox protein Hox-1C
Human Swiss Prot No.:	P20719
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HoxA5 Polyclonal Antibody detects endogenous levels of HoxA5 protein.
Storage and Stability:	-20°C/1 year