

Antibodies Product Datasheet

Product Name:	HoxA6 Antibody
Catalog Number:	C16295
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HOXA6. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	HOXA6
Protein Name:	Homeobox protein Hox-A6
Other Name:	HOXA6; HOX1B; Homeobox protein Hox-A6; Homeobox protein Hox-1B
Human Swiss Prot No.:	P31267
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HoxA6 Polyclonal Antibody detects endogenous levels of HoxA6 protein.
Storage and Stability:	-20°C/1 year