

Antibodies Product Datasheet

Product Name:	HoxB2 Antibody
Catalog Number:	C16300
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HOXB2. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	HOXB2
Protein Name:	Homeobox protein Hox-B2
Other Name:	HOXB2; HOX2H; Homeobox protein Hox-B2; Homeobox protein Hox-2.8; Homeobox protein Hox-2H; K8
Human Swiss Prot No.:	P14652
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HoxB2 Polyclonal Antibody detects endogenous levels of HoxB2 protein.
Storage and Stability:	-20°C/1 year