

Antibodies Product Datasheet

Product Name:	INTS2 Antibody
Catalog Number:	C16387
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human INTS2. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	INTS2
Protein Name:	Integrator complex subunit 2
Other Name:	INTS2; KIAA1287; Integrator complex subunit 2; Int2
Human Swiss Prot No.:	Q9H0H0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	INTS2 Polyclonal Antibody detects endogenous levels of INTS2 protein.
Storage and Stability:	-20°C/1 year