

Antibodies Product Datasheet

Product Name:	ICAM-2 Antibody
Catalog Number:	C16389
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ICAM2. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	ICAM2
Protein Name:	Intercellular adhesion molecule 2
Other Name:	ICAM2; Intercellular adhesion molecule 2; ICAM-2; CD antigen CD102
Human Swiss Prot No.:	P13598
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ICAM-2 Polyclonal Antibody detects endogenous levels of ICAM-2 protein.
Storage and Stability:	-20°C/1 year