

Antibodies Product Datasheet

Product Name:	KLHL29 Antibody
Catalog Number:	C16417
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KLHL29. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	KLHL29
Protein Name:	Kelch-like protein 29
Other Name:	KLHL29; KBTBD9; KIAA1921; Kelch-like protein 29; Kelch repeat and BTB domain-containing protein 9
Human Swiss Prot No.:	Q96CT2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KLHL29 Polyclonal Antibody detects endogenous levels of KLHL29 protein.
Storage and Stability:	-20°C/1 year