

Antibodies Product Datasheet

Product Name:	CD158f2 Antibody
Catalog Number:	C16427
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIR2DL5B. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	KIR2DL5B
Protein Name:	Killer cell immunoglobulin-like receptor 2DL5B
Other Name:	KIR2DL5B; CD158F; CD158F2; KIR2DL5; KIR2DLX; Killer cell immunoglobulin-like receptor 2DL5B; CD158 antigen-like family member F2; Killer cell immunoglobulin-like receptor 2DLX; CD antigen CD158f2
Human Swiss Prot No.:	Q8NHK3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD158f2 Polyclonal Antibody detects endogenous levels of CD158f2 protein.
Storage and Stability:	-20°C/1 year