

Antibodies Product Datasheet

Product Name:	KIF13B Antibody
Catalog Number:	C16437
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIF13B. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	KIF13B
Protein Name:	Kinesin-like protein KIF13B
Other Name:	KIF13B; GAKIN; KIAA0639; Kinesin-like protein KIF13B; Kinesin-like protein GAKIN
Human Swiss Prot No.:	Q9NQT8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KIF13B Polyclonal Antibody detects endogenous levels of KIF13B protein.
Storage and Stability:	-20°C/1 year