

Antibodies Product Datasheet

Product Name:	AASD-PPT Antibody
Catalog Number:	C16447
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AASDHPPT. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	AASDHPPT
Protein Name:	L-aminoadipate-semialdehyde dehydrogenase-phosphopantetheinyl transferase
Other Name:	AASDHPPT; CGI-80; HAH-P; HSPC223; x0005; L-aminoadipate-semialdehyde dehydrogenase-phosphopantetheinyl transferase; 4'-phosphopantetheinyl transferase; Alpha-aminoadipic semialdehyde dehydrogenase-phosphopantetheinyl transferase; AASD-PPT;
Human Swiss Prot No.:	Q9NRN7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AASD-PPT Polyclonal Antibody detects endogenous levels of AASD-PPT protein.
Storage and Stability:	-20°C/1 year