

## Antibodies Product Datasheet

|                               |                                                                                                                            |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | LRRC41 Antibody                                                                                                            |
| <b>Catalog Number:</b>        | C16455                                                                                                                     |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                          |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                     |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/5000. Not yet tested in other applications. |
| <b>Modification:</b>          | Unmodified/Total                                                                                                           |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                          |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.      |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human LRRC41. AA range:391-440                         |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                    |
| <b>Gene Name:</b>             | LRRC41                                                                                                                     |
| <b>Protein Name:</b>          | Leucine-rich repeat-containing protein 41                                                                                  |
| <b>Other Name:</b>            | LRRC41; MUF1; PP7759; Leucine-rich repeat-containing protein 41;<br>Protein Muf1                                           |
| <b>Human Swiss Prot No.:</b>  | Q15345                                                                                                                     |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                    |
| <b>Specificity:</b>           | LRRC41 Polyclonal Antibody detects endogenous levels of LRRC41 protein.                                                    |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                               |