

Antibodies Product Datasheet

Product Name:	ILT-3 Antibody
Catalog Number:	C16465
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LILRB4. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	LILRB4
Protein Name:	Leukocyte immunoglobulin-like receptor subfamily B member 4
Other Name:	LILRB4; ILT3; LIR5; Leukocyte immunoglobulin-like receptor subfamily B member 4; CD85 antigen-like family member K; Immunoglobulin-like transcript 3; ILT-3; Leukocyte immunoglobulin-like receptor 5; LIR-5; Monocyte inhibitory receptor HM18;
Human Swiss Prot No.:	Q8NHJ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ILT-3 Polyclonal Antibody detects endogenous levels of ILT-3 protein.
Storage and Stability:	-20°C/1 year