

Antibodies Product Datasheet

Product Name:	LPD lipase Antibody
Catalog Number:	C16481
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LIPI. AA range:289-338
Concentration:	1 mg/ml
Gene Name:	LIPI
Protein Name:	Lipase member I
Other Name:	LIPI; LPDL; PRED5; Lipase member I; LIPI; Cancer/testis antigen 17; CT17; LPD lipase; Membrane-associated phosphatidic acid-selective phospholipase A1-beta; mPA-PLA1 beta
Human Swiss Prot No.:	Q6XZB0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LPD lipase Polyclonal Antibody detects endogenous levels of LPD lipase protein.
Storage and Stability:	-20°C/1 year