

Antibodies Product Datasheet

Product Name:	ACSVL4 Antibody
Catalog Number:	C16490
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC27A4. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	SLC27A4
Protein Name:	Long-chain fatty acid transport protein 4
Other Name:	SLC27A4; ACSVL4; FATP4; Long-chain fatty acid transport protein 4; FATP-4; Fatty acid transport protein 4; Solute carrier family 27 member 4
Human Swiss Prot No.:	Q6P1M0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ACSVL4 Polyclonal Antibody detects endogenous levels of ACSVL4 protein.
Storage and Stability:	-20°C/1 year