

Antibodies Product Datasheet

Product Name:	LDLRAD3 Antibody
Catalog Number:	C16505
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LDLRAD3. AA range:37-86
Concentration:	1 mg/ml
Gene Name:	LDLRAD3
Protein Name:	Low-density lipoprotein receptor class A domain-containing protein 3
Other Name:	LDLRAD3; Low-density lipoprotein receptor class A domain-containing protein 3
Human Swiss Prot No.:	Q86YD5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LDLRAD3 Polyclonal Antibody detects endogenous levels of LDLRAD3 protein.
Storage and Stability:	-20°C/1 year