

Antibodies Product Datasheet

Product Name:	LRP3 Antibody
Catalog Number:	C16512
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LRP3. AA range:17-66
Concentration:	1 mg/ml
Gene Name:	LRP3
Protein Name:	Low-density lipoprotein receptor-related protein 3
Other Name:	LRP3; Low-density lipoprotein receptor-related protein 3; LRP-3; 105 kDa low-density lipoprotein receptor-related protein; hLRp105
Human Swiss Prot No.:	O75074
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LRP3 Polyclonal Antibody detects endogenous levels of LRP3 protein.
Storage and Stability:	-20°C/1 year