

Antibodies Product Datasheet

Product Name:	LSP1 Antibody
Catalog Number:	C16525
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LSP1. AA range:104-153
Concentration:	1 mg/ml
Gene Name:	LSP1
Protein Name:	Lymphocyte-specific protein 1
Other Name:	LSP1; WP34; Lymphocyte-specific protein 1; 47 kDa actin-binding protein; 52 kDa phosphoprotein; pp52; Lymphocyte-specific antigen WP34
Human Swiss Prot No.:	P33241
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LSP1 Polyclonal Antibody detects endogenous levels of LSP1 protein.
Storage and Stability:	-20°C/1 year