

Antibodies Product Datasheet

Product Name:	DC-LAMP Antibody
Catalog Number:	C16533
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LAMP3. AA range:274-323
Concentration:	1 mg/ml
Gene Name:	LAMP3
Protein Name:	Lysosome-associated membrane glycoprotein 3
Other Name:	LAMP3; DCLAMP; TSC403; Lysosome-associated membrane glycoprotein 3; LAMP-3; Lysosomal-associated membrane protein 3; DC-lysosome-associated membrane glycoprotein; DC LAMP; Protein TSC403; CD antigen CD208
Human Swiss Prot No.:	Q9UQV4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DC-LAMP Polyclonal Antibody detects endogenous levels of DC-LAMP protein.
Storage and Stability:	-20°C/1 year