

Antibodies Product Datasheet

Product Name:	MAD2L1BP Antibody
Catalog Number:	C16545
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MAD2L1BP . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	MAD2L1BP
Protein Name:	MAD2L1-binding protein
Other Name:	MAD2L1BP; CMT2; KIAA0110; MAD2L1-binding protein; Caught by MAD2 protein
Human Swiss Prot No.:	Q15013
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAD2L1BP Polyclonal Antibody detects endogenous levels of MAD2L1BP protein.
Storage and Stability:	-20°C/1 year