

Antibodies Product Datasheet

Product Name:	Med18 Antibody
Catalog Number:	C16578
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED18. AA range:70-119
Concentration:	1 mg/ml
Gene Name:	MED18
Protein Name:	Mediator of RNA polymerase II transcription subunit 18
Other Name:	MED18; Mediator of RNA polymerase II transcription subunit 18; Mediator complex subunit 18; p28b
Human Swiss Prot No.:	Q9BUE0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Med18 Polyclonal Antibody detects endogenous levels of Med18 protein.
Storage and Stability:	-20°C/1 year