

Antibodies Product Datasheet

Product Name:	Rec8 Antibody
Catalog Number:	C16596
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human REC8. AA range:154-203
Concentration:	1 mg/ml
Gene Name:	REC8
Protein Name:	Meiotic recombination protein REC8 homolog
Other Name:	REC8; REC8L1; Meiotic recombination protein REC8 homolog; Cohesin Rec8p
Human Swiss Prot No.:	O95072
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rec8 Polyclonal Antibody detects endogenous levels of Rec8 protein.
Storage and Stability:	-20°C/1 year