

Antibodies Product Datasheet

Product Name:	pICln Antibody
Catalog Number:	C16639
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLNS1A. AA range:184-233
Concentration:	1 mg/ml
Gene Name:	CLNS1A
Protein Name:	Methylosome subunit pICln
Other Name:	CLNS1A; CLCI; ICLN; Methylosome subunit pICln; Chloride channel; nucleotide sensitive 1A; Chloride conductance regulatory protein ICln; I(Cln); Chloride ion current inducer protein; CICI; Reticulocyte pICln
Human Swiss Prot No.:	P54105
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	pICln Polyclonal Antibody detects endogenous levels of pICln protein.
Storage and Stability:	-20°C/1 year