

Antibodies Product Datasheet

Product Name:	MGST2 Antibody
Catalog Number:	C16641
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MGST2 . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	MGST2
Protein Name:	Microsomal glutathione S-transferase 2
Other Name:	MGST2; GST2; Microsomal glutathione S-transferase 2; Microsomal GST-2; Microsomal GST-II
Human Swiss Prot No.:	Q99735
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MGST2 Polyclonal Antibody detects endogenous levels of MGST2 protein.
Storage and Stability:	-20°C/1 year