

Antibodies Product Datasheet

Product Name:	MRP-S25 Antibody
Catalog Number:	C16655
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS25. AA range:124-173
Concentration:	1 mg/ml
Gene Name:	MRPS25
Protein Name:	28S ribosomal protein S25 mitochondrial
Other Name:	MRPS25; RPMS25; 28S ribosomal protein S25; mitochondrial; MRP-S25; S25mt
Human Swiss Prot No.:	P82663
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S25 Polyclonal Antibody detects endogenous levels of MRP-S25 protein.
Storage and Stability:	-20°C/1 year