

Antibodies Product Datasheet

Product Name:	MRP-S27 Antibody
Catalog Number:	C16656
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS27. AA range:16-65
Concentration:	1 mg/ml
Gene Name:	MRPS27
Protein Name:	28S ribosomal protein S27 mitochondrial
Other Name:	MRPS27; KIAA0264; 28S ribosomal protein S27; mitochondrial; MRP-S27; S27mt
Human Swiss Prot No.:	Q92552
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S27 Polyclonal Antibody detects endogenous levels of MRP-S27 protein.
Storage and Stability:	-20°C/1 year