

Antibodies Product Datasheet

Product Name:	MRP-S30 Antibody
Catalog Number:	C16659
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS30. AA range:390-439
Concentration:	1 mg/ml
Gene Name:	MRPS30
Protein Name:	28S ribosomal protein S30 mitochondrial
Other Name:	MRPS30; PDCD9; BM-047; 28S ribosomal protein S30; mitochondrial; MRP-S30; S30mt; Programmed cell death protein 9
Human Swiss Prot No.:	Q9NP92
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S30 Polyclonal Antibody detects endogenous levels of MRP-S30 protein.
Storage and Stability:	-20°C/1 year