

Antibodies Product Datasheet

Product Name:	MRP-S33 Antibody
Catalog Number:	C16661
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPS33. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	MRPS33
Protein Name:	28S ribosomal protein S33 mitochondrial
Other Name:	MRPS33; CGI-139; PTD003; 28S ribosomal protein S33; mitochondrial; MRP-S33; S33mt
Human Swiss Prot No.:	Q9Y291
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-S33 Polyclonal Antibody detects endogenous levels of MRP-S33 protein.
Storage and Stability:	-20°C/1 year