

Antibodies Product Datasheet

Product Name:	MRP-L50 Antibody
Catalog Number:	C16675
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPL50. AA range:109-158
Concentration:	1 mg/ml
Gene Name:	MRPL50
Protein Name:	39S ribosomal protein L50 mitochondrial
Other Name:	MRPL50; 39S ribosomal protein L50; mitochondrial; L50mt; MRP-L50
Human Swiss Prot No.:	Q8N5N7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-L50 Polyclonal Antibody detects endogenous levels of MRP-L50 protein.
Storage and Stability:	-20°C/1 year