

Antibodies Product Datasheet

Product Name:	MPP1 Antibody
Catalog Number:	C16700
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIF20B. AA range:1771-1820
Concentration:	1 mg/ml
Gene Name:	KIF20B
Protein Name:	Kinesin-like protein KIF20B
Other Name:	KIF20B; KRMP1; MPHOSPH1; Kinesin-like protein KIF20B; Cancer/testis antigen 90; CT90; Kinesin-related motor interacting with PIN1; M-phase phosphoprotein 1; MPP1
Human Swiss Prot No.:	Q96Q89
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MPP1 Polyclonal Antibody detects endogenous levels of MPP1 protein.
Storage and Stability:	-20°C/1 year