

Antibodies Product Datasheet

Product Name:	MRP7 Antibody
Catalog Number:	C16727
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRP7. AA range:835-884
Concentration:	1 mg/ml
Gene Name:	ABCC10
Protein Name:	Multidrug resistance-associated protein 7
Other Name:	ABCC10; MRP7; SIMRP7; Multidrug resistance-associated protein 7; ATP-binding cassette sub-family C member 10
Human Swiss Prot No.:	Q5T3U5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP7 Polyclonal Antibody detects endogenous levels of MRP7 protein.
Storage and Stability:	-20°C/1 year