

Antibodies Product Datasheet

Product Name:	ABCC12 Antibody
Catalog Number:	C16728
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRP9. AA range:691-740
Concentration:	1 mg/ml
Gene Name:	ABCC12
Protein Name:	Multidrug resistance-associated protein 9
Other Name:	ABCC12; MRP9; Multidrug resistance-associated protein 9; ATP-binding cassette sub-family C member 12
Human Swiss Prot No.:	Q96J65
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABCC12 Polyclonal Antibody detects endogenous levels of ABCC12 protein.
Storage and Stability:	-20°C/1 year