

Antibodies Product Datasheet

Product Name:	EGFL4 Antibody
Catalog Number:	C16731
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEGF8. AA range:1243-1292
Concentration:	1 mg/ml
Gene Name:	MEGF8
Protein Name:	Multiple epidermal growth factor-like domains protein 8
Other Name:	MEGF8; C19orf49; EGFL4; KIAA0817; Multiple epidermal growth factor-like domains protein 8; Multiple EGF-like domains protein 8; Epidermal growth factor-like protein 4; EGF-like protein 4
Human Swiss Prot No.:	Q7Z7M0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EGFL4 Polyclonal Antibody detects endogenous levels of EGFL4 protein.
Storage and Stability:	-20°C/1 year