

Antibodies Product Datasheet

Product Name:	EGFL5 Antibody
Catalog Number:	C16732
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEGF9. AA range:383-432
Concentration:	1 mg/ml
Gene Name:	MEGF9
Protein Name:	Multiple epidermal growth factor-like domains protein 9
Other Name:	MEGF9; EGFL5; KIAA0818; Multiple epidermal growth factor-like domains protein 9; Multiple EGF-like domains protein 9; Epidermal growth factor-like protein 5; EGF-like protein 5
Human Swiss Prot No.:	Q9H1U4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EGFL5 Polyclonal Antibody detects endogenous levels of EGFL5 protein.
Storage and Stability:	-20°C/1 year