

Antibodies Product Datasheet

Product Name:	MAZ Antibody
Catalog Number:	C16741
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAZ. AA range:159-208
Concentration:	1 mg/ml
Gene Name:	MAZ
Protein Name:	Myc-associated zinc finger protein
Other Name:	MAZ; ZNF801; Myc-associated zinc finger protein; MAZI; Pur-1; Purine-binding transcription factor; Transcription factor Zif87; ZF87; Zinc finger protein 801
Human Swiss Prot No.:	P56270
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAZ Polyclonal Antibody detects endogenous levels of MAZ protein.
Storage and Stability:	-20°C/1 year