

Antibodies Product Datasheet

Product Name:	Myomesin-2 Antibody
Catalog Number:	C16760
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYOM2. AA range:612-661
Concentration:	1 mg/ml
Gene Name:	MYOM2
Protein Name:	Myomesin-2
Other Name:	MYOM2; Myomesin-2; 165 kDa connectin-associated protein; 165 kDa titin-associated protein; M-protein; Myomesin family member 2
Human Swiss Prot No.:	P54296
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Myomesin-2 Polyclonal Antibody detects endogenous levels of Myomesin-2 protein.
Storage and Stability:	-20°C/1 year