

Antibodies Product Datasheet

Product Name:	MYBPC3 Antibody
Catalog Number:	C16776sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MYBPC3 Polyclonal
Concentration:	1 mg/ml
Gene Name:	MYBPC3
Protein Name:	Myosin-binding protein C, cardiac-type (Cardiac MyBP-C) (C-protein, cardiac muscle isoform)
Other Name:	Myosin-binding protein C, cardiac-type (Cardiac MyBP-C) (C-protein, cardiac muscle isoform)
Human Swiss Prot No.:	Q14896
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MYBPC3.
Storage and Stability:	-20°C/1 year